

iTRAQ-based proteomics approaches

Introduced By :


- 1- Habeba Ayman**
- 2- Shorouk Mohamed**
- 3- Eman Mohamed**
- 4- Doaa Tahaa**
- 5- Doaa Abd el-hafez**

❏ Definition:

- **Isobaric tags for relative and absolute quantitation (iTRAQ):** is an isobaric labeling method used in quantitative proteomics by tandem mass spectrometry to determine the amount of proteins from different sources in **a single experiment** .
- It uses stable isotope labeled molecules that can be covalent bonded to the N-terminus and side chain amines of proteins.
- Use multiplexed isobaric mass tags to label peptide digest mixtures.
- Available always are set of **4-plex and 8-plex** mass tags that can be used to label and get quantitative information on up to 4 and 8 different samples in one experiment.
 - Current identification software is ProQuant software.
 - Tracker software has been developed to extract reporter ion peak ratios from non-centroided tandem MS peak lists.
 - Developed by Applied Biosystems in corporation in 2004.


❑ Objective of developing iTRAQ :

- Many differential effects on proteins themselves come from post- translational modifications.
- Protein expression cannot be measured or identified by looking at the mRNA levels.
- By studying effector molecules will contribute to better understanding of disease & in developing new treatment.



❑ The structure of iTRAQ reagents:


- unique charged reporter group.
- Peptide reactive group.
- Neutral balance group.
- The peptide reactive group links an iTRAQ with each lysine side chain and N-terminal group.



❑ The principle of iTRAQ -iTRAQ reagent 4-plex :

- Mixture of peptides is obtained by hydrolyzing the protein sample.
- All peptide in the sample are labeled with different iTRAQ reagents.
- All labeled protein samples are mixed.
- The peptides are subjected to mass spectrometry to obtain the mass spectrum.
- The same peptide from different sources is completely identical in molecular weight and they appear as the same peak samples,


Principle of iTRAQ-iTRAQ reagent 4-plex as an example



This is the first stage of mass spectrometry.


In the second stage of mass spectrometry ,

- The bond between the balance group and peptide reactive group is broken and balance group is lost .
- The same peptide with different isotopic labels produces different masses of reporter ions.
- Reporter ions exhibit different peaks quantitative information of the same peptides between different samples can be obtained by analyzing related data with software.
- Each sample in data bases in the general workflow is reduce system block.
- Each sample digested with trypsin.
- Each sample labeled with a different itraq tag in a single tube and then we can combine all labeled samples into one sample mixture for liquid chromatography tandem mass spectrometry analysis.



❑ Factors affecting results of iTRAQ :

- Evaluation of labeling efficiency and isotope impurity correction .
- Ratio compression and its correction.
- Reporter ion intensity dynamic range.
- Effect of unique and shared peptides in inferring protein ratios.
- Estimation of protein fold changes.



❑ Advantages of iTRAQ :

- Using mass spectrometry .
- Ability to analyze proteins from cell and tissues.
- Multiplexing ability.
- Reduce overall time and variation.
- Cutting-edge facilities & optimized protocols
- Untargeted approach for biomarker discovery

❑ Disadvantages of iTRAQ:

- Costly.
- Very sensitive to contamination from salts.
- A requirement of sophisticated software.
- Variability arising due to the inefficient enzymatic digestion.